

ALCOSCAN

AL2500 SERIES

Highly Sensitive and Reliable Accuracy
Especially useful for self checking the B.A.C.



What difference does AL2500series have?




AL2500
No need to use Mouthpiece




AL2500elite
Use Mouthpiece

Sensor Type



Semiconductor Alcohol sensor

Accuracy




0.40‰
0.50‰
0.60‰

Weight

85g
(Including batteries)




Dimension (mm)



104 (H)
20 (T)
40 (W)

Power



AA size
Alkaline batteries

What Certificate does AL2500 have? CE, DOT, FDA 510K listed